

CHOOSE

PEROXigard[®]

Disinfectants for Life Sciences



POWERFUL

Cleans and disinfects against pathogens of concern in one step

RESPONSIBLE

Non-irritating to eyes and skin, protecting your animals and staff

COMPATIBLE

Non-corrosive formulas, guarding your equipment

FASTER

Effective in as little as 1 minute*, ensuring compliance

PeroxigardDisinfectants.com



Exclusively available at:



PEROXigard®

FAMILY OF PRODUCTS

READY TO USE

Peroxigard® Ready to Use is a unique pre-mixed solution that requires no dilution. Ideal for cleaning and disinfecting biosafety cabinets, laboratory equipment, animal enclosure areas and high touch surfaces.

Formats: 32oz, 1gal, 5gal, 30gal, 55gal

Shelf Life: 2 years (Opened/Unopened)



WIPES

Peroxigard Wipes provide an easy, convenient and effective way to disinfect your surfaces and prevent cross contamination. Ideal for work stations, lab benches, electronic equipment (screens, phones), and other small surfaces.

Formats: 6x7, 160 ct.

Shelf Life: 2 years (Opened/Unopened)

CONCENTRATE

Peroxigard Concentrate is easily diluted and comes in a variety of sizes. A cost effective solution ideal for use on floors, ceilings, walls, cages, kennels and other large areas.

- 1:64 (2oz /gal) Daily Use (Virucidal*)
- 1:16 (8oz /gal) Bactericidal, Virucidal, & Fungicidal

Formats: 1 gal, 5gal, 55gal

Shelf Life: 3 years (Opened/Unopened). Diluted solution remains effective for up to 90 days in a spray bottle.



**See reference sheet for a complete list of efficacy claims, additional features and directions for use.*

IDEAL FOR THE FOLLOWING AREAS

- Animal Examination & Procedure Rooms
- Animal Research Facilities & Vivariums
- Breeding Establishments
- Enclosure & IVC Housing Areas
- Footbaths & Mopping Applications
- Large & Small Animal Facilities
- Laboratories & Equipment
- Non-Human Primate Areas
- Rodent Facilities
- Transportation Vehicles & Equipment